

BCR-ABL1 t(9;22) (mBCR and μ BCR) QUANTITATIVE DETECTION p190/p230

ORDERING INFORMATIONS

REF: *ONC-016-25*
CND Code: *W01060208- t (9;22)*
RDM Code: *1822476/R*
Tests: *25*
Reactions: *50 x 2*
Manufacturer: *BioMol Laboratories s.r.l.*

CONTENTS OF THE KIT

The kit consists of reagents for reverse transcription and PCR amplification.
**the reagents for total RNA extraction are not supplied in the kit*

For in vitro diagnostic use



PRODUCT CHARACTERISTICS

Device belonging to the family of in vitro medical devices REAL-TIME QUANTITATIVE PCR-SOMATIC MUTATIONS. Quantitative detection of BCR-ABL1 t(9;22) breakpoint m-bcr (e1a3 e e1a2) and μ -bcr (e18a2, e18a3, e19a2 e e19a3) transcripts by reverse transcription, amplification with oligonucleotides and specific probes and subsequent detection with qPCR-Real-time using plasmids for standard curve. Kit optimized for Real-Time PCR instrumentation Biorad CFX96 Dx, Biorad Opus Dx, Agilent AriaDx.

SCIENTIFIC BACKGROUND

Myeloproliferative neoplasms (MPNs) are hematologic malignancies characterized by the proliferation of one or more myeloid lineages: granulocytic, erythroid, megakaryocytic, and/or mast cell. According to the 2016 World Health Organization criteria, the MPN classification includes seven subcategories: chronic myeloid leukemia (CML), chronic neutrophilic leukemia, polycythemia vera (PV), primary myelofibrosis (PMF), essential thrombocythemia (ET), eosinophilic leukemia chronic - not otherwise specified and MPN, unclassifiable (MPN-U).

The Philadelphia chromosome (Ph) derived from the translocation between chromosomes 9 and 22 with subsequent BCR-ABL1 fusion, is present in about 95% of cases of chronic myeloid leukemia (CML), in 25-30% of cases of acute lymphoblastic leukemia (ALL) of adults and in 2-4% of ALL of children.

§ Genetic basis and molecular pathophysiology of classical myeloproliferative neoplasms. *Blood*. 2017 Feb 9; 129(6):667-679. doi: 10.1182/blood-2016-10-695940. Epub 2016 Dec 27. Review.

§ Classification and diagnosis of myeloproliferative neoplasms: the 2008 World Health Organization criteria and point-of-care diagnostic algorithms. *Leukemia*. 2008 Jan; 22(1):14-22. Epub 2007 Sep 20. Review.

§ The 2016 revision to the World Health Organization classification of myeloid neoplasms and acute leukemia. *Blood*. 2016 May 19; 127(20): 2391-405. Epub 2016 Apr 11.

§ Guidelines for the measurement of BCR-ABL1 transcripts in chronic myeloid leukaemia. *Br J Haematol*. 2011 Apr; 153(2):179-90. doi: 10.1111/j.1365-2141.2011.08603.x. Epub 2011 Mar 8.

§ Leukemia. 2003 Dec;17(12):2318-57. doi: 10.1038/sj.leu.2403135. Standardization and quality control studies of 'real-time' quantitative reverse transcriptase polymerase chain reaction of fusion gene transcripts for residual disease detection in leukemia - a Europe Against Cancer program.

§ J Clin Oncol. 2009 Dec 10;27(35):6041-51. doi: 10.1200/JCO.2009.25.0779. Epub 2009 Nov 2. Chronic myeloid leukemia: an update of concepts and management recommendations of European LeukemiaNet


§ Leukemia. 2009 Nov;23(11):1957-63. doi: 10.1038/leu.2009.168. Epub 2009 Aug 27. Harmonization of molecular monitoring of CML therapy in Europe

§ European LeukemiaNet (2009). Chronic myeloid leukemia: an update of concepts and management recommendations of European LeukemiaNet. *Journal of Clinical Oncology*, 27, 6041-6051.

§ Leukemia. 2015 May;29(5):999-1003. doi: 10.1038/leu.2015.29. Epub 2015 Feb 5. Laboratory recommendations for scoring deep molecular responses following treatment for chronic myeloid leukemia

CLINICAL SIGNIFICANCE

The BCR-ABL1 rearrangement results in the generation of fusion proteins with constitutive tyrosine kinase activity. Based on the specific breakpoints of the rearrangement, different isoforms of the BCR-ABL1 fusion protein are generated, which correlate with different leukemic phenotypes. Three breakpoint regions in the BCR gene have been described: major (M-BCR), minor (m-BCR), and micro (μ -BCR). More than 95% of Ph+ CML patients have the rearrangement in the M-BCR region (p210 BCR-ABL1), with the e13a2 and e14a2 transcripts most represented. The breakpoint in the m-BCR region generates the p190 BCR-ABL1 protein with the e1a2 transcript mostly represented. A third BCR-ABL1 protein, p230 BCR-ABL1 (μ BCR), can also be observed. This translocation is associated with CML characterized by granulocytic hyperplasia and, in general, with a more indolent clinical course.



BCR-ABL1 t(9;22) (mBCR e μ BCR) QUANTITATIVE DETECTION p190/p230

ORDERING INFORMATIONS

REF: *ONC-016-25*
 CND Code: *W01060208- t (9;22)*
 RDM Code: *1822476/R*
 Tests: *25 x2*
 Reactions: *100*
 Manufacturer: *BioMol Laboratories s.r.l.*

CONTENTS OF THE KIT

The kit consists of: reagents for reverse transcription and PCR amplification.
**the reagents for total RNA extraction are not supplied in the kit*

For in vitro diagnostic use



CONTENTS OF THE KIT

| DESCRIPTION | LABEL | LABEL | VOLUME | STORAGE |
|---|---|---|-------------------|---------|
| | | | ONC-016-25 | |
| Mix oligonucleotides and probes | Mix PCR p190 BCR-ABL1 4X | | 1 x 250 μ l | - 20 °C |
| Mix oligonucleotides and probes | Mix PCR p230 BCR-ABL1 4X | | 1 x 250 μ l | - 20 °C |
| Mix buffer and RT/Taq polymerase enzyme | Mix RT-PCR 4X | | 1 x 500 μ l | - 20 °C |
| Deionized H ₂ O | Deionized H ₂ O | | 1 x 1 ml | - 20 °C |
| Recombinant DNA | CAL 1 p190/abl - 1,08 x10 ⁵ copies | CAL 1 p230/abl - 1,08 x10 ⁵ copies | 1 x 30 μ l | - 20 °C |
| Recombinant DNA | CAL 2 p190/abl -1,08 x10 ⁴ copies | CAL 2 p230/abl -1,08 x10 ⁴ copies | 1 x 30 μ l | - 20 °C |
| Recombinant DNA | CAL 3 p190/abl -1,08 x10 ³ copies | CAL 3 p230/abl -1,08 x10 ³ copies | 1 x 30 μ l | - 20 °C |
| Recombinant DNA | CAL 4 p190/abl - 1,08 x10 ² copies | CAL 4 p230/abl - 1,08 x10 ² copies | 1 x 30 μ l | - 20 °C |
| Recombinant DNA | CAL 5 p190/abl - 10,8 copies | CAL 5 p230/abl - 10,8 copies | 1 x 30 μ l | - 20 °C |
| Recombinant RNA | Positive control p190/p230/abl | | 1 x 60 μ l | - 20 °C |
| Recombinant RNA | Negative control housekeeping | | 1 x 60 μ l | - 20 °C |

TECHNICAL CHARACTERISTICS

COD. ONC-016-25

| | |
|---|---|
| STABILITY | 18 months |
| REAGENTS STATUS | Ready to use |
| BIOLOGICAL MATRIX | Total RNA extracted from white blood cells from whole blood or bone marrow aspirate |
| POSITIVE AND NEGATIVE CONTROLS | Recombinant RNA for 3 analytical sessions |
| STANDARD CURVE | Recombinant DNA p190 and p230, 5 points with known concentration from 10,8 to 1,085 copies, plasmid standard curve |
| TECHNOLOGY | RT-PCR ONE STEP in Real-time; oligonucleotides and specific probes; 2 FAM/HEX fluorescence channels |
| VALIDATED INSTRUMENTS | Biorad CFX96 Dx, Biorad Opus Dx and Agilent AriaDx |
| RUNNING TIME | 110 min |
| THERMAL CYCLING PROFILE | 1 cycle at 25 °C (2 min); 1 cycle at 50 °C (25 min); 1 cycle at 95 °C (2 min); 50 cycles at 95 °C (5 sec) + 60 °C (45 sec). Reading at 60°C |
| ANALYTICAL SPECIFICITY | Absence of non-specific pairings of oligonucleotides and probes; absence of cross-reactivity |
| LIMIT OF DETECTION (LOD) | = 10 copies |
| LIMIT OF BLANK (LOB) | 0% NCN |
| REPRODUCIBILITY | 99,9% |
| DIAGNOSTIC SPECIFICITY / DIAGNOSTIC SENSITIVITY | 100%/98% |